



ROTAID™
AED CABINETS



ROTAID LIST
AED Overview



Serial number **HEATING PLATE**


The serial number of the Rotaid AED cabinet can be found here:





PLUS & SOLID PLUS (HEAT LED) AED PLACEMENT ASSISTANT

Use both bolts to hang the AED. Please check your brand and type to find the ideal position.



**2x DIN screw
M6 required**

Cardiac Science
Powerheart G3 Plus
Powerheart G3 Pro
Powerheart G5

Defibtech
Lifeline
Lifeline Pro
Lifeline View

HeartSine
Samaritan Pad 350P
Samaritan Pad 360P
Samaritan Pad 450P

Mediana
A10
A15

Mindray
BeneHeart D1
BeneHeart D1 Pro

Nihon Kohden
Cardiolife AED-2150
Cardiolife AED-3100

Philips
Heartstart FR3
Heartstart FRx
Heartstart OnSite/
HS1

Physio control
CR2
Lifepak 1000
Lifepak CR Plus
Lifepak Express

Schiller
Fred Easy
Fred Easy Life
Fred Easyport
Fred PA-1

**2x DIN screw
M6 required**

CU Medical Systems
I-Pad CU-SP1
I-Pad NF1200

Telefunken
Telefunken AED

**Special holding
bracket required
(install without bag)**

Zoll
AED 3
AED Plus

**2x DIN screw
M6 required**

Primedic
Heartsave AED
Heartsave AS

**Special bracket
required**


Cardiaid
CT0207RF/S

Defibtech
Lifeline View

Primedic
Heartsave One
Heartsave Pad

Progetti
Rescue Sam
Rescue Sam Pro

2x DIN screw M6
Required for: — — — —





ROTAID AED LIST

AED List

* Click on your AED and go directly to the correct page: 



Cardiac Science

Powerheart G3 Plus
Powerheart G3 Pro
Powerheart G5



Philips

Heartstart FR3
Heartstart FRx
Heartstart OnSite/HS1



Cardiaid

CT0207RF
CT0207RS



Physio control

CR2
Lifepak 1000
Lifepak CR Plus
Lifepak Express



Medical Systems, Inc.

CU Medical Systems

I-Pad CU-SP1
I-Pad NF1200



Defibtech

Lifeline
Lifeline Pro
Lifeline View



Primedic

Heartsave AED
Heartsave AS
Heartsave One
Heartsave Pad



HeartSine

Samaritan Pad 350P
Samaritan Pad 360P
Samaritan Pad 450P



Progetti

Rescue Sam
Rescue Sam Pro



Mediana

A10
A15



Schiller

Fred Easy
Fred Easy Life
Fred Easyport
Fred PA-1



Mindray

BeneHeart D1
BeneHeart D1 Pro



Telefunken

Telefunken AED



Nihon Kohden

Cardiolife AED-2150
Cardiolife AED-3100



Zoll (install without bag)

AED 3
AED Plus



Cardiac Science
Powerheart G3





Cardiaid
CT0207R (F/S)





CU Medical Systems
I-Pad (CU-SP1/NF1200)





Defibtech **Lifeline** (Pro)





Defibtech Lifeline View





HeartSine
Samaritan PAD (350P/360P/450P)





Mediana
A10





Mediana
A15





Mindray
Beneheart D1 (Pro)





Nihon Kohden
CardioLife AED-2100





Nihon Kohden
CardioLife AED-3100





Philips
Heartstart (HS1/FRx)



HS1



FRx



Physio Control CR2





Physio Control
Lifepak (1000/CR Plus/Express)





Pimedec
HeartSave (Aed/AS)





Pimedec
Heartsave (One/Pad)





Progetti
Rescue Sam (PRO)





Schiller
Fred Easy (Life/-port)





Telefunken
Telefunken





Zoll
AED 3
(install without bag)





Zoll
AED Plus
(install without bag)

